

# User Manual

## **1. Position Yourself in Front of the Camera**

Sit directly in front of the webcam so that your face and eyes are clearly visible. Make sure the room has good lighting for accurate eye detection.

## **2. Start the Application**

Run the program on the computer. Three windows will appear: Frame Window (webcam feed), Virtual Keyboard Window (keyboard letters), and Board Window (typed text).

## **3. Select the Keyboard Side**

At the start, two options will appear: LEFT and RIGHT. Look left to select the left keyboard and look right to select the right keyboard. Hold your gaze briefly to confirm the selection.

## **4. Keyboard Layout**

The keyboard is divided into two sections: Left Side and Right Side.

| Left Side Keyboard | Right Side Keyboard |
|--------------------|---------------------|
| Q W E R T          | Y U I O P           |
| A S D F G          | H J K L             |
| Z X C V            | B N M               |
| < (Backspace)      | _ (Space)           |
|                    | ENT (New Line)      |

## **5. Letter Selection**

After selecting a keyboard side, letters will highlight one by one automatically. Wait until the desired letter is highlighted.

## **6. Blink to Type**

When the required letter is highlighted, blink your eyes to select it. The selected letter will appear on the Board window.

## **7. Special Keys**

\_ : Space, < : Backspace (delete previous character), ENT : New line.

### **8. *Continue Typing***

After typing a letter, the system will again ask you to choose LEFT or RIGHT keyboard. Repeat the process to continue typing.

### **9. *Exit the System***

Press the ESC key on the keyboard to close the application.